

“

NEWMALLE

”

Anatomic Ankle Locking Plate & Screw System

Fibular Locking Plate, Lateral (Lt, Rt)

3.5mm Locking Screw, Gold

Fibular Locking Plate, Posterior (Lt, Rt)

3.5mm Cortical Screw

Distal Tibia Locking Plate, Anterior (Lt, Rt)

Distal Tibia Locking Plate, Posterior (Lt, Rt)

Distal Tibia Locking Plate, Medial

Distal Tibia Locking Plate, Medial, Hook

Variable Locking Calcaneus Plate (50, 60, 70mm)

2.7mm Locking Screw, Lt Blue

2.7mm Cortical Screw, Low Profile Head



TDM - TRAUMA

NEWMALLE



Contents

NEWMALLE

- TDM Principles
- TDM Fracture Solution
- Developing Background
- KEY Features
- Surgical Technique
- Ordering Information Part I.
- Ordering Information Part II.

TDM PRINCIPLES



It has already been more than a decade since the company was founded to become a leading company for fracture treatment. While keeping the basics, all employees have worked together and have been quick to reflect customers' needs.

The constant dedication and efforts of physicians for patients always present us with new challenges.

TDM knows that these challenges are none other than the minds of doctors for patients.

To know this, we wanted to reflect customer requests into our products and services as much as possible, and we will continue to do so in the future.

“

TDM Principle is simple.


Basic Plus One.

That **ONE** for the **patients.**

”

TDM Fracture Treatment Solution

An increase in osteoporosis patients in the elderly and a fracture caused by sports activities in the younger population means a rise in fractures with bone defects.



Implant

We keep pace with the growing number of elderly patients and the trend of trauma from sports activities. Based on the Basic plus one principle.

Biologics

**Fracture
Healing
System**

We approach patients systematically for their quick return to daily life.

Variable Locking Foot Plate & Screw Developing Background

Surgeon's Demand

Design for accurate reconstruction and keep the anatomic shape for the healing period.



“

NEWMALLE

”

- **EXTREME THIN**
- **ANATOMICALLY SHAPED**
- **SYNDESMOSIS HOLES**
- **LOCKING PLATE SYSTEM**


For the patients

Minimize discomfort of patients during treatment.

TDM R&D

It will be as thin and structurally rigid as possible, while at the same time being able to make various attempts.

KEY Features:



Product Features

- Anatomical Fit
- Variety of designs according to the bone anatomy and fracture type: 6 type plates
- Optimal fracture fixation of Trimalleolar , Tibia And Fibular Fractures
- Syndesmosis Holes repairs a high ankle sprain and provides stabilization of a fibular fractures.
- Locking and Combi hole options
- Low-profile Plates creates minimum irritation
- Applicable for Osteotomies and Nonunions



Surgical Technique Guide



TDM as the manufacturer of this device, does not practice medicine and does not recommend this or any other surgical technique for use on a specific patient. The surgeon who performs any implant procedure is responsible for determining and using the appropriate techniques for implanting the device in each patient.

Surgical Tech. Newmalle – Fibular Locking Plate

01

Preparation

Place the patient in a supine position with an ipsilateral pelvic sandbag or in a mid-lateral position.

Using a lateral approach, make a straight incision to enable the best view of the lateral talar and calcaneal walls, mobilizing the peroneal tendons. (Figure 1)

An appropriate size of the plate is selected from the system.

As the plate is positioned, place the plate on the bone directly or fixate the plate position with a Guide Pin. (Figure 2)



Figure 1




Figure 2

Surgical Tech. Newmalle – Fibular Locking Plate

02

Preparation

Attach the Drilling Block onto the plate. (Figure3)

Drill the screw hole through $\text{\O}2.0$ Locking Drill Sleeve by using $\text{\O}2.0\text{mm}$ Drill bit. (Figure 4)

Note:

An alternative option for $\text{\O}2.0\text{mm}$ Locking Drill Sleeve $\text{\O}2.0\text{mm}$ is $\text{\O}2.0\text{mm}$ Drill Sleeve Single handle.

Instruments

Drill Bit

901-01220 $\text{\O} 2.0\text{mm} \times 130\text{mm}$

Drilling Block, Fibular, Lateral

Left 901-20268

Right 901-20269

Locking Drill Sleeve, Single Hand

901-02220, $\text{\O}2.0\text{mm}$

Locking Drill Sleeve

901-12520, $\text{\O}2.0 \times 60\text{mm}$





Figure 3



Figure 4

Surgical Tech. – Fibular Locking Plate

03

Placement


Measure and determine the screw length for the distal and the shaft parts with a depth gauge. (Figure 5)

Instrument

Depth Gauge 2.7/3.5/4.0
901-04027



Figure 5



Surgical Tech. – Fibular Locking Plate

04

Depth and Screw insertion

Select and pick up a corresponding screw from the measurement.

Insert the screw by using a corresponding size of the screwdriver.

-2.7mm Locking Screw for the Distal part.


-3.5mm Screw for the Shaft part.

(Figure 7, 8)


Note:

[Checking Reduction and Fixation]

An intraoperative radiograph is recommended to check the position of the screws and the final reduction of the fracture.



2.7 Locking Screw



3.5 Locking Screw

Instrument

Screwdriver, Holding

901-35020, Hex. 2.0mm

901-35025, Hex. 2.5mm



Figure 7




Figure 8

Ordering Information, Implant


NEWMALLE Plate, Fibular

Fibular Locking Plate, Lateral

166-13103 Lt/3H
166-13104 Lt/4H
166-13105 Lt/5H
166-13107 Lt/7H

166-13203 Rt/3H
166-13204 Rt/4H
166-13205 Rt/5H
166-13207 Rt/7H

2.7mm / 3.5mm Screw (Hole)




Fibular Locking Plate, Lateral

Fibular Locking Plate, Posterior

166-14103 Lt/3H
166-14104 Lt/4H
166-14105 Lt/5H

166-14203 Lt/3H
166-14204 Lt/4H
166-14205 Lt/5H

2.7mm / 3.5mm Screw (Hole)



Fibular Locking Plate, Posterior

Ordering Information, Implant


NEWMALLE Plate, Distal Tibia

Distal Tibia Locking Plate, Anterior

160-11103 Lt/3H
160-11104 Lt/4H
160-11105 Lt/5H
160-11106 Lt/6H

160-11203 Rt/3H
160-11204 Rt/4H
160-11205 Rt/5H
160-11206 Rt/6H

2.7mm/ 3.5mm Screw




Hole

Distal Tibia Locking Plate, Posterior

161-14103 Lt/3H
161-14104 Lt/4H
161-14105 Lt/5H
161-14106 Lt/6H

161-14203 Rt/3H
161-14204 Rt/4H
161-14205 Rt/5H
161-14206 Rt/6H


2.7mm / 3.5mm Screw (Hole)



Distal Tibia Locking Plate, Medial & Medial Hook


161-15003 3H
161-15004 4H

3.5mm Screw



161-16003 3H
161-16004 4H

3.5mm Screw



Hole

Ordering Information, Implant


NEWMALLE Plate, Calcaneus

Variable Locking Calcaneus

144-11050 50mm

144-11060 60mm

144-11070 70mm





?

Ordering Information, Implant

2.7mm Screw (2.0 Hex)

2.7mm Locking Screw, Lt Blue

227-26010 10mm
227-26012 12mm
227-26014 14mm
227-26016 16mm
227-26018 18mm
227-26020 20mm
227-26022 22mm
227-26024 24mm
227-26026 26mm
227-26028 28mm
227-26030 30mm
227-26032 32mm
227-26034 34mm
227-26036 36mm
227-26038 38mm
227-26040 40mm



2.7mm Cortical Screw, Low Profile Head

227-27010 10mm
227-27012 12mm
227-27014 14mm
227-27016 16mm
227-27018 18mm
227-27020 20mm
227-27022 22mm
227-27024 24mm
227-27026 26mm
227-27028 28mm
227-27030 30mm
227-27032 32mm
227-27034 34mm
227-27036 36mm
227-27038 38mm
227-27040 40mm

Ordering Information, Implant

3.5mm Screw (2.5 Hex)

3.5mm Locking Screw, Gold

- 235-26010 10mm
- 235-26012 12mm
- 235-26014 14mm
- 235-26016 16mm
- 235-26018 18mm
- 235-26020 20mm
- 235-26022 22mm
- 235-26024 24mm
- 235-26026 26mm
- 235-26028 28mm
- 235-26030 30mm
- 235-26032 32mm
- 235-26034 34mm
- 235-26036 36mm
- 235-26038 38mm
- 235-26040 40mm
- 235-26045 45mm **Option**
- 235-26050 50mm **Option**



3.5mm Cortical Screw

- 235-21010 10mm
- 235-21012 12mm
- 235-21014 14mm
- 235-21016 16mm
- 235-21018 18mm
- 235-21020 20mm
- 235-21022 22mm
- 235-21024 24mm
- 235-21026 26mm
- 235-21028 28mm
- 235-21030 30mm
- 235-21032 32mm
- 235-21034 34mm
- 235-21036 36mm
- 235-21038 38mm
- 235-21040 40mm
- 235-21045 45mm **Option**
- 235-21050 50mm **Option**
- 235-21055 55mm **Option**
- 235-21060 60mm **Option**

Ordering Information, Instrument

Drill Bit

901-01220, Ø 2.0 x 130mm

901-01125, Ø 2.5 x 130mm / Option

901-22127, Ø 2.7 x 130mm

901-22135, Ø 3.5 x 160mm / Option




Locking Drill Sleeve

901-12520, Ø2.0x60mm

901-12527, Ø2.7x60mm

901-12220, 2.0x60mm / Option

901-12227, 2.7x60mm / Option



Locking Drill Sleeve, Single Handle

901-02220, Ø2.0mm

901-02227, Ø2.7mm

901-13325, Ø2.5mm / Option

901-13335, Ø3.5mm / Option




Depth Gauge

901-04027 2.7/3.5/4.0mm



Holding Forcep / Option

901-08001



Drill Bit, AO Chuck / Option

901-31120 Ø2.0x 130mm

901-31025, Ø2.5x130mm

901-31127, Ø2.7x130mm

901-31135, Ø2.7x130mm



Screwdriver, Holding

901-35020, Hex. 2.0mm

901-35025, Hex. 2.5mm



Screwdriver Shaft, Holding

901-15420 Hex 2.0mm/Small Chuck

901-15425 Hex 2.5mm/Small Chuck



Torque Limit Handle

901-17008, 0.8Nm / Small Chuck

901-17015, 1.5Nm / Small Chuck



Bender

901-06035, Left

901-06036, Right



Ordering Information, Instrument

Drilling Block, Fibular, Lateral, Left
901-20268

Drilling Block, Fibular, Posterior, Left
901-21268

Drilling Block, Distal Tibia, Anterior, Left
901-22268

Drilling Block, Distal Tibia, Lateral, Left
901-23268




Drilling Block, Fibular, Lateral, Right
901-20269

Drilling Block, Fibular, Posterior, Right
901-21269

Drilling Block, Distal Tibia, Anterior, Right
901-22269

Drilling Block, Distal Tibia, Lateral, Right
901-23269

Ordering Information, Tray



Container for Ankle Locking Implants

Outer Cover 999103-0001S

Outer Body 999103-0002S

Plate Rack 999103-0003S

* In this container, there are plates & corresponding Drilling Blocks

Trifix Universal Screw & Instruments

Universal Container

Outer Cover 901-99301S

Outer Body 901-99302S

Outer Cover 901-99306S


Outer Cover 901-99307A

Outer Cover 901-99308S

* In this container, there are screws & corresponding instruments.



Ordering Information, Tray



Container for Ankle Locking Implants

Outer Cover 999103-0001S

Outer Body 999103-0002S

Plate Rack 999103-0003S

* In this container, there are plates & corresponding Drilling Blocks

Trifix Universal Screw & Instruments

Universal Container

Outer Cover 901-99301S

Outer Body 901-99302S

Outer Cover 901-99306S

Outer Cover 901-99307A

Outer Cover 901-99308S

* In this container, there are screws & corresponding instruments.



Ordering Information, Other Products



1.0mm Mini Cable System

It is designed to be easy to use but also to function when necessary. It can be applied to a variety of areas. **Hand, Foot, Clavicle, Fibular, Ulna, Patella...** Not only is it used alone, but it can also be used in conjunction with a variety of plates to produce the optimum effect.



Headless Compression Cannulated System (HCC)

Designed to minimize soft tissue irritation through Headless Fixation
Compression achieved along the length of the screw due to its continuously
Self-Cutting , Self-tapping, Helical Relief Flutes
The opposing surface thread design optimizes compression toward the fracture/osteotomy
Available in 2.3/3.5, 3.0/4.0, 3.5/5.0mm screw diameters



Threaded Lock Pin System, K-Wire, Steinman Pin and Suture Wire

TDM is providing basic and diverse products for fracture treatment.

TDM - TRAUMA

NEWMALLE

“NEWMALLE”

Follow us On



www.tradimedics.com/



www.facebook.com/tdmkorea/



www.youtube.com/tdmkorea



www.instagram.com/tdmcoltd/



linkedin.com/company/42851518/



Head Office

69, Cheomdan Venture so-ro
37beon-gil, Buk-gu, Gwangju, Korea (61003)
Tel. +82-31-732-0632 Fax. +82-70-4941-3236